



The
right
connection
is **everything.**

Medical Fittings

Product Catalog



Medical Fittings Product Catalog



Overview.....	4
Ohmeda Style Quick Connect.....	6
Ohmeda Style Winged Quick Connect.....	7
Compact Ohmeda Style Check Unit (Female Coupler).....	8
Ohmeda Style Check Unit (Female Coupler).....	9
Chemetron Style Quick Connect.....	10
Chemetron Style Check Unit (Female Coupler)	11
DISS Female.....	12
DISS Male.....	14
Puritan-Bennett Style Quick Connect	16
Puritan-Bennett Style Check Unit (Female Coupler)	17
Schrader Style Quick Connect	18
Schrader Style Check Unit (Female Coupler).....	19
Medstar Style Quick Connect.....	20
Oxequip Style Quick Connect	21
“Y” Block Assemblies and DISS “Y” Assemblies	22
“T” Bar Assemblies and DISS “T” Assemblies	23
T-Branch, Y-Block, Hex/Quad Block and Rail Blocks.....	26
Extension Manifolds, Pole Clamp and Ferrules.....	27
Medical Grade Hose	28
Tubing Nipples	29
NPT Pipe Fittings	30
Condensation Traps.....	32
Neb Drives and Rescue Flows.....	33
International Fittings.....	34



The right connection is everything.


Quality products to make the right connection.

Precision Medical is a global pioneer in respiratory care focused on developing and manufacturing high-quality medical devices. With three main product lines spanning the hospital and homecare markets, our product portfolio covers infants through advanced age and ranges from high-quality oxygen administration to aspiration and suction products.

Precision Medical's durable medical fittings, which are made in the USA, allow full customization and targeting to meet the needs of our customers and users. Our fittings are found in hospitals, surgery centers, and long-term/sub-acute facilities worldwide. We manufacture adapters to fit every need such as both male and female couplers and quick connects for all major styles including Ohmeda, Chemetron, Puritan-Bennett, Schrader, and DISS connections. Customized configurations and hose assemblies are available to meet your unique requirements.

At Precision Medical, we take great pride in providing exceptional customer service. Our knowledgeable representatives are here to help you match the right product to your specific need and have it shipped quickly. We are 100% committed to providing products that are safe, accurate, durable and easy to maintain.

Precision Medical. The Difference is Precision.



Ohmeda Style Quick Connect



Oxygen



Air















Vacuum



Nitrous Oxide



WAGD

	Oxygen	Air	Vacuum	Nitrous Oxide	WAGD
 1/8" NPT Male	2101	2201	2301	2401	2601
 1/8" NPT Male x 1.5" Long			2301X		
 1/8" NPT Female	2102	2202	2302	2402	2602
 1/4" NPT Male	2103	2203	2303	2403	2603
 DISS Male (O2)	2105				
 DISS Male		2205	2305	2405	2605
 DISS Hex Nut	2106	2206	2306	2406	2606
 DISS Hand Tight	2108	2208	2308	2408	2608
 Hose Barb	1/4" Hose Barb 2111	1/4" Hose Barb 2211	1/4" Hose Barb 2311 5/16" Hose Barb 2312	1/4" Hose Barb 2411	1/4" Hose Barb 2611 5/16" Hose Barb 2615
 Hose Barb Swivel	1/4" Hose Barb 2111S	1/4" Hose Barb 2211S	1/4" Hose Barb 2311S 5/16" Hose Barb 2312S	1/4" Hose Barb 2411S	1/4" Hose Barb 2611S 5/16" Hose Barb 2615S
 DISS Male with Check Valve (O2)	2114				
 DISS Male with Check Valve		2214	2314	2414	2614

Ohmeda Style Winged Quick Connect



Oxygen



Air



Vacuum



Nitrous Oxide

	Oxygen	Air	Vacuum	Nitrous Oxide
 1/8" NPT Male	2101W	2201W	2301W	2401W
 1/8" NPT Male x 1.5" Long			2301WX	
 1/8" NPT Female	2102W	2202W	2302W	2402W
 1/4" NPT Male	2103W	2203W	2303W	2403W
 DISS Male (O2)	2105W			
 DISS Male		2205W	2305W	2405W
 DISS Hex Nut	2106W	2206W	2306W	2406W
 DISS Hand Tight	2108W	2208W	2308W	2408W
 Hose Barb	1/4" Hose Barb 2111W	1/4" Hose Barb 2211W	1/4" Hose Barb 2311W 5/16" Hose Barb 2312W	1/4" Hose Barb 2411W
 Hose Barb Swivel	1/4" Hose Barb 2111WS	1/4" Hose Barb 2211WS	1/4" Hose Barb 2311WS 5/16" Hose Barb 2312WS	1/4" Hose Barb 2411WS
 DISS Male with Check Valve (O2)	2114W			
 DISS Male with Check Valve		2214W	2314W	2414W

Compact Ohmeda Style Check Unit

(Female Coupler)



Oxygen



Air



Vacuum



Nitrous Oxide



WAGD



1/8" NPT Female

2160

2260

2360

2460

2660



1/8" NPT Male

2161

2261

2361

2461

2661



1/4" NPT Male

2163

2263

2363

2463

2663



Hose Barb

1/4" Hose Barb
2164

1/4" Hose Barb
2264

1/4" Hose Barb
2364
5/16" Hose Barb
2367

1/4" Hose Barb
2464

1/4" Hose Barb
2664
5/16" Hose Barb
2667



DISS Male (O2)

2165



DISS Male

2265

2365

2465

2665



DISS Hex Nut

2166

2266

2366

2466

2666



DISS Hand Tight

2168

2268

2368

2468

2668

Ohmeda Style Check Unit

(Female Coupler)

	Oxygen	Air	Vacuum	Nitrous Oxide	WAGD
 1/8" NPT Female	2150	2250	2350	2450	2650
 1/8" NPT Male	2151	2251	2351	2451	2651
 1/4" NPT Male	2153	2253	2353	2453	2653
 Hose Barb	1/4" Hose Barb 2154	1/4" Hose Barb 2254	1/4" Hose Barb 2354 5/16" Hose Barb 2357	1/4" Hose Barb 2454	1/4" Hose Barb 2654 5/16" Hose Barb 2657
 DISS Male (O2)	2155				
 DISS Male		2255	2355	2455	2655
 DISS Hex Nut	2156	2256	2356	2456	2656
 DISS Hand Tight	2158	2258	2358	2458	2658

Chemetron Style Quick Connect



Oxygen



Air















Vacuum



Nitrous Oxide



WAGD

	Oxygen	Air	Vacuum	Nitrous Oxide	WAGD
 1/8" NPT Male	3101	3201	3301	3401	3601
 1/8" NPT Male x 1.5" Long			3301X		
 1/8" NPT Female	3102	3202	3302	3402	3602
 1/4" NPT Male	3103	3203	3303	3403	3603
 DISS Male (O2)	3105				
 DISS Male		3205	3305	3405	3605
 DISS Hex Nut	3106	3206	3306	3406	3606
 DISS Hand Tight	3108	3208	3308	3408	3608
 Hose Barb	1/4" Hose Barb 3111	1/4" Hose Barb 3211	1/4" Hose Barb 3311 5/16" Hose Barb 3312	1/4" Hose Barb 3411	1/4" Hose Barb 3611 5/16" Hose Barb 3615
 Hose Barb Swivel	1/4" Hose Barb 3111S	1/4" Hose Barb 3211S	1/4" Hose Barb 3311S 5/16" Hose Barb 3312S	1/4" Hose Barb 3411S	1/4" Hose Barb 3611S 5/16" Hose Barb 3615S
 DISS Male with Check Valve (O2)	3114				
 DISS Male with Check Valve		3214	3314	3414	3614

Chemetron Style Check Unit

(Female Coupler)



Oxygen



Air



Vacuum



Nitrous Oxide









WAGD

	Oxygen	Air	Vacuum	Nitrous Oxide	WAGD
	3150	3250	3350	3450	3650
1/8" NPT Female					
	3151	3251	3351	3451	3651
1/8" NPT Male					
	3153	3253	3353	3453	3653
1/4" NPT Male					
	1/4" Hose Barb 3154	1/4" Hose Barb 3254	1/4" Hose Barb 3354 5/16" Hose Barb 3357	1/4" Hose Barb 3454	1/4" Hose Barb 3654 5/16" Hose Barb 3657
Hose Barb					
	3155				
DISS Male (O2)					
		3255	3355	3455	3655
DISS Male					
	3156	3256	3356	3456	3656
DISS Hex Nut					
	3158	3258	3358	3458	3658
DISS Hand Tight					













DISS Female

	Oxygen	Air	Vacuum	Nitrous Oxide	Nitrogen	WAGD	Carbon Dioxide
 Hex Nut-1/8" NPT with 1" Stem	0106						
 Hex Nut-1/8" NPT with 1" Stem Swivel	0106S						
 Hex Nut-1/8" NPT with 1/2" Stem	0107	0207	0307	0407	0507		504248
 Hex Nut-1/4" NPT with 1/2" Stem					0506	0606	
 Hex Nut-Hose Barb	1/4" Hose Barb 0115	1/4" Hose Barb 0216	1/4" Hose Barb 0316 5/16" Hose Barb 0317	1/4" Hose Barb 0416	1/4" Hose Barb 0516	1/4" Hose Barb 0616 5/16" Hose Barb 0617	1/4" Hose Barb 504154
 Hex Nut-Hose Barb Swivel with 1/2" Stem	1/4" Hose Barb 0115S	1/4" Hose Barb 0216S	1/4" Hose Barb 0316S 5/16" Hose Barb 0317S	1/4" Hose Barb 0416S	1/4" Hose Barb 0516S	1/4" Hose Barb 0616S 5/16" Hose Barb 0617S	1/4" Hose Barb 504154S
 Oxygen DISS Female x Oxygen DISS Female	0104						

DISS Female

	Oxygen	Air	Vacuum	Nitrous Oxide	Nitrogen	WAGD	Carbon Dioxide
 Hand Tight- 1/8" NPT with 1" Stem	0110						
 Hand Tight- Hose Barb Swivel	0110S						
 Hand Tight- 1/8" NPT with 1/2" Stem	0111	0209	0309	0409	0509		507883
 Hand Tight- 1/4" NPT with 1/2" Stem					0508	0608	
 Hand Tight- Hose Barb	1/4" Hose Barb 0118	1/4" Hose Barb 0218	1/4" Hose Barb 0318 5/16" Hose Barb 0319	1/4" Hose Barb 0418	1/4" Hose Barb 0518	1/4" Hose Barb 0618 5/16" Hose Barb 0619	1/4" Hose Barb 504287
 Hand Tight- Hose Barb Swivel with 1/2" Stem	1/4" Hose Barb 0118S	1/4" Hose Barb 0218S	1/4" Hose Barb 0318S 5/16" Hose Barb 0319S	1/4" Hose Barb 0418S	1/4" Hose Barb 0518S	1/4" Hose Barb 0618S 5/16" Hose Barb 0619S	1/4" Hose Barb 504287S

DISS Male

							
	Oxygen	Air	Vacuum	Nitrous Oxide	Nitrogen	WAGD	Carbon Oxide
DISS Male to 1/8" NPT Male	0701	0731	0751	0771	0791		
DISS Male to 1/4" NPT Male	0703	0733	0753	0773	0792	0794	8946
 1/4" Hose Barb	0709					0601 Demand Check 0602	
 1/4" Hose Barb with Knob	0710	0740	0760	0780	0795		
 1/4" Hose Barb with Knob & Check	0711	0741	0761	0781	0796		
 1/4" Hose Barb with Knob Swivel	0710S	0740S	0760S	0780S	0795S		
 1/4" Hose Barb with Knob & Check Swivel	0711S	0741S	0761S	0781S	0796S		



Oxygen



Air



Vacuum



Nitrous Oxide



Nitrogen



WAGD



Carbon Oxide

Demand Check Valve

DISS Male with Demand Check Valve to 1/8" NPT Male	0712	0742		0782	0797		
DISS Male with Demand Check Valve to 1/4" NPT Male	0713	0743	0763	0783	0798	0799	8947

One Way Check Valve

	O2, DISS Male x 1/4" NPT Male	508373
	Air, DISS Male x 1/4" NPT Male	508374
	1/4" NPT Female x 1/4" NPT Male	508375

<p>0705 0706</p> <p>Oxygen DISS 90° Elbow Male x Hex Nut</p>	<p>1" Stem 0705</p> <p>1½" Stem 0706</p>	<p>Oxygen Power Take-Off (1/8" NPT)</p>	0150
<p>Oxygen DISS Male x DISS Male</p>	0704	<p>Air Power Take-Off (1/8" NPT)</p>	0151
<p>Oxygen DISS Male x 1/8" NPT Female</p>	0702		

Puritan-Bennett Style Quick Connect



Oxygen



Air















Vacuum



Nitrous Oxide



WAGD

	Oxygen	Air	Vacuum	Nitrous Oxide	WAGD
 1/8" NPT Male	4101	4201	4301	4401	4601
 1/8" NPT Male x 1.5" Long			4301X		
 1/8" NPT Female	4102	4202	4302	4402	4602
 1/4" NPT Male	4103	4203	4303	4403	4603
 DISS Male with Knob (O2)	4105				
 DISS Male with Knob		4205	4305	4405	4605 (no knob)
 DISS Hex Nut	4106	4206	4306	4406	4606
 DISS Hand Tight	4108	4208	4308	4408	4608
 Hose Barb with Knob	1/4" Hose Barb 4111	1/4" Hose Barb 4211	1/4" Hose Barb 4311 5/16" Hose Barb 4312 (no knob)	1/4" Hose Barb 4411	1/4" Hose Barb 4611 (no knob) 5/16" Hose Barb 4615 (no knob)
 Hose Barb Swivel	1/4" Hose Barb 4111S	1/4" Hose Barb 4211S	1/4" Hose Barb 4311S 5/16" Hose Barb 4312S	1/4" Hose Barb 4411S	1/4" Hose Barb 4611S 5/16" Hose Barb 4615S
 DISS Male with Knob and Check Valve (O2)	4114				
 DISS Male with Knob and Check Valve		4214	4314	4414	4614 (no knob)

Puritan-Bennett Style Check Unit

(Female Coupler)



Oxygen



Air



Vacuum

	Oxygen	Air	Vacuum
	4150	4250	4350
	4151	4251	4351
	4153	4253	4353
	1/4" Hose Barb 4154	1/4" Hose Barb 4254	1/4" Hose Barb 4354 5/16" Hose Barb 4357
	4155		
		4255	4355
	4156	4256	4356
	4158	4258	4358

Schrader Style Quick Connect



Oxygen

Air

Vacuum

Nitrous Oxide



1/8" NPT Male

6101

6201

6301

6401



1/8" NPT Male x 1.5" Long

6301X



1/8" NPT Female

6102

6202

6302

6402



1/4" NPT Male

6103

6203

6303

6403



DISS Male (O2)

6105



DISS Male

6205

6305

6405



DISS Hex Nut

6106

6206

6306

6406



DISS Hand Tight

6108

6208

6308

6408



Hose Barb

1/4" Hose Barb
6111

1/4" Hose Barb
6211

1/4" Hose Barb
6311

1/4" Hose Barb
6411



DISS Male
and Check Valve (O2)

6114



DISS Male
and Check Valve

6214

6314

6414

Schrader Piston
Identification



Hex Piston
(standard)



Swivel Piston
(add "S" after part number)

Schrader Style Check Unit

(Female Coupler)



Oxygen



Air



Vacuum



Nitrous Oxide



Air/Nitrogen



1/8" NPT Female

6150

6250

6350

6450




1/8" NPT Male

6151

6251

6351

6451



1/4" NPT Male

6153

6253

6353

6453



Hose Barb

1/4" Hose Barb
6154

1/4" Hose Barb
6254

1/4" Hose Barb
6354

1/4" Hose Barb
6454



DISS Male (O2)

6155



DISS Male

6255

6355

6455



DISS Hex Nut

6156

6256

6356

6456

Schrader Air
Coupler x
Nitrogen DISS

6556



DISS Hand Tight

6158

6258

6358

6458

Schrader Air
Coupler x
Nitrogen DISS

6558

Medstar Style Quick Connect



Oxygen



Air



Vacuum

	Oxygen	Air	Vacuum
 1/8" NPT Male	8101	8201	8301
 1/8" NPT Male x 1.5" Long			8301X
 1/8" NPT Female	8102	8202	8302
 1/4" NPT Male	8103	8203	8303
 DISS Male (O2)	8105		
 DISS Male		8205	8305
 DISS Hex Nut	8106	8206	8306
 DISS Hand Tight	8108	8208	8308
 Hose Barb	1/4" Hose Barb 8111	1/4" Hose Barb 8211	1/4" Hose Barb 8311
 Hose Barb Swivel	8111S	8211S	8311S
 DISS Male with Check Valve (O2)	8114MS		
 DISS Male with Check Valve		8214	8314

Oxequip Style Quick Connect

	Oxygen	Air	Vacuum
 1/8" NPT Male	5101	5201	5301
 1/8" NPT Male x 1.5" Long			5301X
 1/8" NPT Female	5102	5202	5302
 1/4" NPT Male	5103	5203	5303
 DISS Male with Knob	5105		
 DISS Male with Knob		5205	5305
 DISS Hex Nut	5106	5206	5306
 DISS Hand Tight	5108	5208	5308
 Hose Barb	1/4" Hose Barb 5111	1/4" Hose Barb 5211	1/4" Hose Barb 5311
 DISS Male with Knob and Check Valve	5114		
 DISS Male with Knob and Check Valve		5214	5314

“Y” Block Assemblies and DISS “Y” Assemblies



Oxygen



Air






Vacuum



Nitrous Oxide



Nitrogen

	Oxygen	Air	Vacuum	Nitrous Oxide	Nitrogen
 Ohmeda Style	2152Y	2252Y	2352Y	2452Y	
with Winged Quick Connect	2152WY	2252WY	2352WY	2352WY	
 Compact Ohmeda Style	2162Y	2262Y	2362Y	2462Y	
with Winged Quick Connect	2162WY	2262WY	2362WY	2462WT	
 Chemetron Style	3152Y	3252Y	3352Y	3452Y	
 Puritan-Bennett Style	4152Y	4252Y	4352Y		
 Schrader Style	6152Y	6252Y	6352Y	6452Y	
 DISS “Y” Assemblies	7211	7212	7213	7214	7240

“T” Bar Assemblies and DISS “T” Assemblies

	 Oxygen	 Air	 Vacuum	 Nitrous Oxide
 Ohmeda Style	2152T	2252T	2352T	2452T
with Winged Quick Connect	2152WT	2252WT	2352WT	2452WT
 Ohmeda Style Offset Right	2152T-R	2252T-R	2352T-R	2452T-R
with Winged Quick Connect	2152WT-R	2252WT-R	2352WT-R	2452WT-R
 Ohmeda Style Offset Left	2152T-L	2252T-L	2352T-L	2452T-L
with Winged Quick Connect	2152WT-L	2252WT-L	2352WT-L	2452WT-L
 Compact Ohmeda Style	2162T	2262T	2362T	2462T
with Winged Quick Connect	2162WT	2262WT	2362WT	2462WT
 Compact Ohmeda Style Offset Right	2162T-R	2262T-R	2362T-R	2462T-R
with Winged Quick Connect	2162WT-R	2262WT-R	2362WT-R	2462WT-R
 Compact Ohmeda Style Offset Left	2162T-L	2262T-L	2362T-L	2462T-L
with Winged Quick Connect	2162WT-L	2262WT-L	2362WT-L	2462WT-L

“T” Bar Assemblies and DISS “T” Assemblies

	 Oxygen	 Air	 Vacuum	 Nitrous Oxide
 Chemetron Style	3152T	3252T	3352T	3452T
 Chemetron Style Offset Right	3152T-R	3252T-R	3352T-R	3452T-R
 Chemetron Style Offset Left	3152T-L	3252T-L	3352T-L	3452T-L
 Puritan-Bennett Style	4152T	4252T	4352T	
 Puritan-Bennett Style Offset Right	4152T-R	4252T-R	4352T-R	
 Puritan-Bennett Style Offset Left	4152T-L	4252T-L	4352T-L	

“T” Bar Assemblies and DISS “T” Assemblies

	 Oxygen	 Air	 Vacuum	 Nitrous Oxide
 <p>Schrader Style</p>	6152T	6252T	6352T	6452T
 <p>Schrader Style Offset Right</p>	6152T-R	6252T-R	6352T-R	6452T-R
 <p>Schrader Style Offset Left</p>	6152T-L	6252T-L	6352T-L	6452T-L
 <p>DISS “T” Assemblies</p>	7111T	7112T		7114T
 <p>DISS “T” Assemblies Offset Right</p>	7111T-R	7112T-R		7114T-R
 <p>DISS “T” Assemblies Offset Left</p>	7111T-L	7112T-L		7114T-L

T-Branch, Y-Block, Hex/Quad Block and Rail Blocks



T-Branch



Y-Block



Hex Block (3 outlet)



Quad Block (4 outlet)

1/8" NPT Standard
7100

1/8" NPT Offset
7100L/R

1/8" NPT Female
7201

1/4" NPT Female
7200

1/8" NPT Female
9054

1/4" NPT Female
9052

1/4" NPT Female
9051

Oxygen

Air

Vacuum



Hex Block
(1/8" NPT Female Inlet)

7217

7218

7219



Quad Block Assembly
(1/8" NPT Female Inlet)

7220

7221

7222


Rail Block Assemblies

All Rail Block Assemblies include
Wall Bracket and Plate

7223

7224

7225



Rail Block Assembly
(DISS Male Inlet)

Replacement Parts

- 7226 — Mounting Bracket
- 9056 — Mounting Plate
- 9054 — Hex Block, 3 outlets
- 9058 — Screw (2 required)

Custom assemblies available upon request.

Extension Manifolds, Pole Clamp and Ferrules

Extension Manifold

FEATURES:

Economical expansion of single wall outlet.

Permanently mounts to wall.

2 foot conductive hose with 1/8" NPT Female inlet is standard.


Outlets are DISS Male with Check Valves.



Oxygen



Air



Vacuum

	Oxygen	Air	Vacuum
 Extension Manifold-2 Outlet	7230	7233	
 Extension Manifold-2 Outlet with Swivel			7236
 Extension Manifold-3 Outlet	7231	7234	
 Extension Manifold-4 Outlet	7232	7235	
 Oxygen Manifold	505568		

Pole Clamp

FEATURES:

1/8" NPT Female holes (top and front) - Comes with 1 plug on bottom



7227A shown



7227B shown

Pole Clamp

1" Pole
7216

1-1/2" Pole
7227

Ferrule



Ferrule

1/4" Hose		5/16" Hose	
3/4" Long	1" Long	3/4" Long	1" Long
9030	507712	503780	507501

Add -25 or -50 at end of Part Number for those quantities









Medical Grade Hose

FEATURES:









Kink resistant and light weight.
 FDA medical grade conductive hose.
 Reinforced with tough nylon braiding.
 Conveniently packaged in 250 foot rolls or as noted.

Color coded to conform to CGA-C9 standards.
 Maximum working pressure to 200 psi @ 70°F.
 Meets NFPA 99 requirements for conductivity (conductive hose only).
 Keep out of sunlight and in a climate controlled atmosphere

Non-Conductive Hose

Oxygen	Green		(1/4" ID-250' Roll)	9020-250
Air	Yellow		(1/4" ID-250' Roll)	9021-250
Vacuum	White		(1/4" ID-250' Roll)	9022-250
Nitrous Oxide	Blue		(1/4" ID-250' Roll)	9023-250
Nitrogen	Black		(1/4" ID-250' Roll)	9028-250
Carbon Dioxide	Gray		(1/4" ID-250' Roll)	9033-250
WAGD	Purple		(5/16" ID-250' Roll)	9038-250
	Clear		(1/4" ID-100' Roll)	504202-100

Conductive Hose

Oxygen	Green		(1/4" ID-250' Roll)	9120-250
Air	Yellow		(1/4" ID-250' Roll)	9121-250
Vacuum	White		(1/4" ID-250' Roll) (5/16" ID-250' Roll)	9122-250 9135-250
Nitrous Oxide	Blue		(1/4" ID-250' Roll)	9123-250
Nitrogen	Black		(1/4" ID-250' Roll)	9128-250
Carbon Dioxide	Gray		(1/4" ID-250' Roll)	9125-250
WAGD	Purple		(1/4" ID-250' Roll)	9124-250
ISO Air	Black & White Striped		(1/4" ID-250' Roll)	9126-250

Custom hose assembly configurator is available at: <http://hoseconfig.precisionmedical.com>



Crimping Tool



Hose Retractor








9130

1/4" Single Clamp	HR-1
1/4" Double Clamp	HR-2
5/16" Single Clamp	HR-3
5/16" Double Clamp	HR-4

Tubing Nipples

Part Number	Method of Cleaning								
	Autoclave	Ethylene Oxide	Chlorine Dioxide	Hydrogen Peroxide	Gamma Ray	Electron Beam	Vaporized Hydrogen Peroxide (VHP)	Quaternary Cleaners	Sodium Hypochlorite "Bleach Solutions"
0154/0154C		X	X	X	X	X	X	X	X
0158/0158C		X	X	X	X	X	X	X	X
0159/0159C	X	X	X	X	X	X	X	X	X
0159K/0159KC	X	X	X	X	X	X	X	X	X
0115	X	X	X	X	X	X	X	X	X








* Please Remove O-Ring for Autoclaving *
 * Clean all Tubing Nipples after each patient use *

 <p>0154 0154C</p>	<p>Oxygen</p> <p>Oxygen with Chain</p>	<p>0154</p> <p>0154C</p>
<p>Hand Tight Tubing Nipple</p>		
 <p>0158 0158C</p>	<p>Oxygen</p> <p>Oxygen with Chain</p>	<p>0158</p> <p>0158C</p>
<p>Wing Nut Tubing Nipple</p>		
 <p>0159 0159C</p>	<p>Oxygen</p> <p>Oxygen with chain</p>	<p>0159</p> <p>0159C</p>
<p>Hex Nut Tubing Nipple</p>		
 <p>0159K 0159KC</p>	<p>Oxygen</p> <p>Oxygen with chain</p>	<p>0159K</p> <p>0159KC</p>
<p>Knurled Nut Tubing Nipple</p>		
	<p>Clear</p> <p>Yellow</p> <p>Black</p> <p>White</p> <p>Green</p>	<p>0119</p> <p>0120</p> <p>0121</p> <p>0122</p> <p>0124</p>
<p>Disposable Tubing Nipple</p>		<p>*Note: sold in packs of 25 and 50. Add -25 or -50 at end of part number for those quantities</p>
 <p>9073 9090</p>	<p>Straight</p> <p>90°</p>	<p>9073</p> <p>9090</p>
<p>Tubing Nipple (for Vacuum Regulator)</p>		
	<p>6" Chain with Closed Ring</p> <p>6" Chain with Open Ring</p>	<p>8004</p> <p>8005</p>
<p>Chain Ring Assembly</p>		

NPT Pipe Fittings

	1/8" NPT	1/4" NPT		1/8" NPT	1/4" NPT
 Pipe Tee	7030	7031	 Plug	7013	7012
 Street Tee	7032	7033	 45° Street Elbow	7056	7057
 Female Tee	7034	7035	 90° Elbow	7058	7059
 Plug	7019	7017	 90° Street Elbow	7063	7065
 Drop Ear	Tee 7081	Elbow 7082	 Master Flow Valve	7064	7045
 Elbow	7092		 Cross	7036	7037

NPT Pipe Fittings

 <p>Hex Nipple</p>	<p>1/8" NPT Male x 1/8" NPT Male 1/4" NPT Male x 1/4" NPT Male 1/8" NPT Male x 1/4" NPT Male</p>	<p>7020 7021 7022</p>
 <p>Connector</p>	<p>1/8" NPT Female x 1/8" NPT Male 1/8" NPT Female x 1/4" NPT Male 1/4" NPT Female x 1/8" NPT Male 1/4" NPT Female x 1/4" NPT Male 1/8" NPT Female x 1/8" NPT Male x 2.5" Long</p>	<p>1017 7024 7025 7026 8955</p>
 <p>Coupler</p>	<p>1/8" NPT Female x 1/8" NPT Female 1/4" NPT Female x 1/4" NPT Female 1/4" NPT Female x 1/8" NPT Female</p>	<p>7027 7028 7029</p>
 <p>Hose Barb Connector</p>	<p>1/8" NPT Male x 1/8" Diameter 1/8" NPT Male x 3/16" Diameter 1/8" NPT Male x 1/4" Diameter 1/4" NPT Male x 1/4" Diameter 1/8" NPT Male x 5/16" Diameter 1/8" NPT Female x 1/4" Diameter 1/4" NPT Female x 1/4" Diameter 1/4" NPT Female x 5/16" Diameter</p>	<p>7038 7039 7040 7043 8952 7046 7049 8951</p>
 <p>Swiveling Hose Barbs</p>	<p>1/8" NPT Male x 5/16" Diameter 1/8" NPT Male x 1/4" Diameter</p>	<p>507765 507764</p>
 <p>Pipe Nipple</p>	<p>1/8" NPT Male x 0.75" Long 1/8" NPT Male x 1.5" Long 1/8" NPT Male x 2" Long 1/8" NPT Male x 2.5" Long 1/4" NPT Male x 1.5" Long 1/4" NPT Male x 2" Long 1/4" NPT Male x 2.5" Long 1/8" NPT Male x 3" Long</p>	<p>7044 7050 7051 7052 7053 7054 7055 7062</p>
 <p>Cap</p>	<p>1/8" NPT Female</p>	<p>7018</p>
 <p>Hose Splicer</p>	<p>1/4" Inside Diameter</p>	<p>7041</p>

Condensation Traps








	Inlet DISS Female Oxygen x Outlet DISS Male Oxygen	PM15-4
	Inlet DISS Female Air x Outlet DISS Male Air	PM15-41
	Inlet DISS Female Air x Outlet DISS Male Air 90° Down	PM15-43
	Inlet DISS Female Oxygen 90° Up x Outlet DISS Male Oxygen 90° Down	PM15-44
	Inlet DISS Male Air 90° Down x Outlet DISS Female Air 90° Up	PM15-45

Replacement Parts



	507588
Jar	
	507587
Filter	
	1271
Jar and Filter	

Neb Drives and Rescue Flows








Neb Drives

		Oxygen	Air
 Ohmeda	6 to 7 l/min	8126	8226
	7 to 8 l/min	8127	8227
	8 to 9 l/min	8128	8228
	10 to 11 l/min	8129	8229
 Chemetron	6 to 7 l/min	8136	8236
	7 to 8 l/min	8137	8237
	8 to 9 l/min	8138	8238
	10 to 11 l/min	8139	8239
 DISS Hand Tight	6 to 7 l/min	8116	8216
	7 to 8 l/min	8117	8217
	8 to 9 l/min	8118	8218
	10 to 11 l/min	8119	8219
 Puritan-Bennett with Knob	6 to 7 l/min	8146	8246
	7 to 8 l/min	8147	8247
	8 to 9 l/min	8148	8248
	10 to 11 l/min	8149	8249
 Schrader	6 to 7 l/min	8166	8266
	7 to 8 l/min	8167	8267
	8 to 9 l/min	8168	8268
	10 to 11 l/min	8169	8269
 Oxequip with Knob	6 to 7 l/min	8156	8256
	7 to 8 l/min	8157	8257
	8 to 9 l/min	8158	8258
	10 to 11 l/min	8159	8259
 Medstar	6 to 7 l/min	8176	8276
	7 to 8 l/min	8177	8277
	8 to 9 l/min	8178	8278
	10 to 11 l/min	8179	8279

Rescue Flows

	Oxygen
	14 to 16 l/min
	DISS Hand Tight 8115
	Ohmeda 8125
	Chemetron 8135
	Puritan-Bennett 8145
	Oxequip 8155
	Schrader 8165
Medstar 8175	
	Oxygen
	24 to 26 l/min
	DISS Hand Tight 8114
	Ohmeda 8124
	Chemetron 8134
	Puritan-Bennett 8144
	Oxequip 8154
	Schrader 8164
Medstar 8174	

International Fittings

	Oxygen	Air	Vacuum	Nitrous Oxide
FRENCH (AFNOR)				
 Probe x 1/8" NPT Male	505979	505980	505981	
 Probe x 1/4" Hose Barb	505766	505767	505768	
EUROPEAN (NIST)				
 Nut/Nipple x 1/8" NPT Male	505752	505772	505773	
 Nut/Nipple x 1/4" Hose Barb	505774	505775	505776	
 Body Adaptor x 1/8" NPT Male	505777	505778	505779	
 Body Adaptor x 1/4" Hose Barb	505780	505781	505782	
AUSTRALIAN (SIS)				
 Nut/Nipple x 1/8" NPT Male	505769	505770	505771	

International Fittings

	Oxygen	Air	Vacuum	Nitrous Oxide
SCANDINAVIAN (AGA)				
 Probe x 1/8" NPT Male	505763	505764	505765	
 Probe x Hose Barb Swivel	1/4" Hose Barb 505763S	1/4" Hose Barb 505764S	1/4" Hose Barb 505765S	
GERMAN (DIN)				
 Probe x 1/8" NPT Male	505753	505761	505762	
 Probe x Hose Barb Swivel	1/4" Hose Barb 505753S	1/4" Hose Barb 505761S	1/4" Hose Barb 505762S	
BRITISH (BOC)				
 Probe x 1/8" NPT Male	6180	6184	6182	6183
 Probe x Hose Barb Swivel	1/4" Hose Barb 6180S	1/4" Hose Barb 6184S	1/4" Hose Barb 6182S	1/4" Hose Barb 6183S
 PK Schrader	6130	6230		



PrecisionMedical

The Difference is Precision.

300 Held Drive | Northampton, PA 18067

precisionmedical.com

T: 800-272-7285 | F: 800-352-1240