

# EasyPulse POC

OXYGEN CONCENTRATOR

The EasyPulse POC package includes:



5-year warranty on concentrator and sieve beds.

External battery and backpack are also available for purchase.



Talk to your provider about insurance coverage or purchase options for you EasyPulse POC!

**PrecisionMedical**

300 Held Drive • Northampton, PA 18067  
800-272-7285 • precisionmedical.com



508644 Rev 3 8/20 Printed in USA

**EasyPulse POC**  
OXYGEN CONCENTRATOR



precisionmedical.com



## Independence and greater mobility

EasyPulse Portable Oxygen Concentrator was carefully designed to give you independence and greater mobility. Whether shopping for groceries, visiting friends and family or going on vacation, EasyPulse POC helps you get back the freedom to go where you want, when you want. Based on your doctor's prescribed oxygen flow, Precision Medical offers EasyPulse in both 3 and 5 liter configurations.

## The highest oxygen purity

Only Precision Medical provides vacuum pressure swing adsorption technology (VPSA) ensuring each sieve bed is perfectly rejuvenated every cycle. This results in clean sieve beds, so they can produce high oxygen purity for a long period of time.

## Tested for purity

Each EasyPulse POC goes through rigorous testing to ensure oxygen purity and durability. Our quality control ensures that every EasyPulse POC leaves the factory producing the highest oxygen purity possible.

## Fixed flow control

EasyPulse POCs are built with advanced Controlled Minute Volume (CMV) technology. No matter your breathing rate, Precision Medical's proprietary CMV technology consistently delivers a fixed amount of oxygen every minute during the first third of inspiration.



Makes traveling with oxygen easier

No more dragging tanks or change overs

Designed to deliver the highest therapeutic oxygen

One touch start and stop

3 or 5 pulse flow setting options

Battery options to increase time away from home

Available with a carry bag or backpack for hands free wear



FAA approved for air travel



Great for cruise ships

## Battery Duration

at any breath rate

EasyPulse POC-3 (PM4130)		EasyPulse POC-5 (PM4150)	
Setting 1	5.5 hours	Setting 1	4.7 hours
Setting 2	3.8 hours	Setting 2	3.4 hours
Setting 3	3.0 hours	Setting 3	2.7 hours
		Setting 4	2.0 hours
		Setting 5	1.6 hours

Battery durations can be doubled by using the optional external battery.