



The Difference is **Precision.**



MEDICAL FITTINGS

Ohmeda Style Quick Connect



Oxygen



Air



Vacuum



Nitrous Oxide



WAGD

	Oxygen	Air	Vacuum	Nitrous Oxide	WAGD
	2101	2201	2301	2401	2601
1/8" NPT Male					
			2301X		
1/8" NPT Male x 1.5" Long					
	2102	2202	2302	2402	2602
1/8" NPT Female					
	2103	2203	2303	2403	2603
1/4" NPT Male					
	2105				
DISS Male (O2)					
		2205	2305	2405	2605
DISS Male					
	2106	2206	2306	2406	2606
DISS Hex Nut					
	2108	2208	2308	2408	2608
DISS Hand Tight					
	1/4" Hose Barb 2111	1/4" Hose Barb 2211	1/4" Hose Barb 2311 5/16" Hose Barb 2312	1/4" Hose Barb 2411	1/4" Hose Barb 2611 5/16" Hose Barb 2615
Hose Barb					
	1/4" Hose Barb 2111S	1/4" Hose Barb 2211S	1/4" Hose Barb 2311S 5/16" Hose Barb 2312S	1/4" Hose Barb 2411S	1/4" Hose Barb 2611S 5/16" Hose Barb 2615S
Hose Barb Swivel					
	2114				
DISS Male with Check Valve (O2)					
		2214	2314	2414	2614
DISS Male with Check Valve					

MEDICAL FITTINGS

Ohmeda Style Winged Quick Connect



Oxygen



Air



Vacuum



Nitrous Oxide

	Oxygen	Air	Vacuum	Nitrous Oxide
	2101W	2201W	2301W	2401W
			2301WX	
	2102W	2202W	2302W	2402W
	2103W	2203W	2303W	2403W
	2105W			
		2205W	2305W	2405W
	2106W	2206W	2306W	2406W
	2108W	2208W	2308W	2408W
	1/4" Hose Barb 2111W	1/4" Hose Barb 2211W	1/4" Hose Barb 2311W 5/16" Hose Barb 2312W	1/4" Hose Barb 2411W
	1/4" Hose Barb 2111WS	1/4" Hose Barb 2211WS	1/4" Hose Barb 2311WS 5/16" Hose Barb 2312WS	1/4" Hose Barb 2411WS
	2114W			
		2214W	2314W	2414W

Compact Ohmeda Style Check Unit (Female Coupler)



Oxygen



Air



Vacuum



Nitrous Oxide



WAGD

	Oxygen	Air	Vacuum	Nitrous Oxide	WAGD
	2160	2260	2360	2460	2660
	2161	2261	2361	2461	2661
	2163	2263	2363	2463	2663
	1/4" Hose Barb 2164	1/4" Hose Barb 2264	1/4" Hose Barb 2364 5/16" Hose Barb 2367	1/4" Hose Barb 2464	1/4" Hose Barb 2664 5/16" Hose Barb 2667
	2165				
		2265	2365	2465	2665
	2166	2266	2366	2466	2666
	2168	2268	2368	2468	2668

MEDICAL FITTINGS

Ohmeda Style Check Unit (Female Coupler)



Oxygen



Air



Vacuum



Nitrous Oxide



WAGD

	Oxygen	Air	Vacuum	Nitrous Oxide	WAGD
 1/8" NPT Female	2150	2250	2350	2450	2650
 1/8" NPT Male	2151	2251	2351	2451	2651
 1/4" NPT Male	2153	2253	2353	2453	2653
 Hose Barb	1/4" Hose Barb 2154	1/4" Hose Barb 2254	1/4" Hose Barb 2354 5/16" Hose Barb 2357	1/4" Hose Barb 2454	1/4" Hose Barb 2654 5/16" Hose Barb 2657
 DISS Male (O2)	2155				
 DISS Male		2255	2355	2455	2655
 DISS Hex Nut	2156	2256	2356	2456	2656
 DISS Hand Tight	2158	2258	2358	2458	2658

Chemetron Style Quick Connect



Oxygen



Air



Vacuum



Nitrous Oxide



WAGD

MEDICAL FITTINGS

	Oxygen	Air	Vacuum	Nitrous Oxide	WAGD
	3101	3201	3301	3401	3601
1/8" NPT Male					
			3301X		
1/8" NPT Male x 1.5" Long					
	3102	3202	3302	3402	3602
1/8" NPT Female					
	3103	3203	3303	3403	3603
1/4" NPT Male					
	3105				
DISS Male (O2)					
		3205	3305	3405	3605
DISS Male					
	3106	3206	3306	3406	3606
DISS Hex Nut					
	3108	3208	3308	3408	3608
DISS Hand Tight					
	1/4" Hose Barb 3111	1/4" Hose Barb 3211	1/4" Hose Barb 3311 5/16" Hose Barb 3312	1/4" Hose Barb 3411	1/4" Hose Barb 3611 5/16" Hose Barb 3615
Hose Barb					
	1/4" Hose Barb 3111S	1/4" Hose Barb 3211S	1/4" Hose Barb 3311S 5/16" Hose Barb 3312S	1/4" Hose Barb 3411S	1/4" Hose Barb 3611S 5/16" Hose Barb 3615S
Hose Barb Swivel					
	3114				
DISS Male with Check Valve (O2)					
		3214	3314	3414	3614
DISS Male with Check Valve					

Chemetron Style Check Unit (Female Coupler)



Oxygen



Air



Vacuum



Nitrous Oxide



WAGD

	3150	3250	3350	3450	3650
1/8" NPT Female					
	3151	3251	3351	3451	3651
1/8" NPT Male					
	3153	3253	3353	3453	3653
1/4" NPT Male					
	1/4" Hose Barb 3154	1/4" Hose Barb 3254	1/4" Hose Barb 3354 5/16" Hose Barb 3357	1/4" Hose Barb 3454	1/4" Hose Barb 3654 5/16" Hose Barb 3657
Hose Barb					
	3155				
DISS Male (O2)					
		3255	3355	3455	3655
DISS Male					
	3156	3256	3356	3456	3656
DISS Hex Nut					
	3158	3258	3358	3458	3658
DISS Hand Tight					

DISS Female Hex Nut

OXYGEN

AIR

VACUUM

NITROUS
OXIDE

NITROGEN

WAGD

CARBON
DIOXIDE

MEDICAL FITTINGS

	0106						
	0106S						
	0107	0207	0307	0407	0507		504248
					0506	0606	
	1/4" Hose Barb 0115	1/4" Hose Barb 0216	1/4" Hose Barb 0316 5/16" Hose Barb 0317	1/4" Hose Barb 0416	1/4" Hose Barb 0516	1/4" Hose Barb 0616 5/16" Hose Barb 0617	1/4" Hose Barb 504154
	1/4" Hose Barb 0115S	1/4" Hose Barb 0216S	1/4" Hose Barb 0316S 5/16" Hose Barb 0317S	1/4" Hose Barb 0416S	1/4" Hose Barb 0516S	1/4" Hose Barb 0616S 5/16" Hose Barb 0617S	1/4" Hose Barb 504154S
	0104						

DISS Female Hand Tight

	OXYGEN	AIR	VACUUM	NITROUS OXIDE	NITROGEN	WAGD	CARBON DIOXIDE
 Hand Tight- 1/8" NPT with 1" Stem	0110						
 Hand Tight- Hose Barb Swivel	0110S						
 Hand Tight- 1/8" NPT with 1 1/2" Stem	0111	0209	0309	0409	0509		507883
 Hand Tight- 1/4" NPT with 1 1/2" Stem					0508	0608	
 Hand Tight- Hose Barb	1/4" Hose Barb 0118	1/4" Hose Barb 0218	1/4" Hose Barb 0318 5/16" Hose Barb 0319	1/4" Hose Barb 0418	1/4" Hose Barb 0518	1/4" Hose Barb 0618 5/16" Hose Barb 0619	1/4" Hose Barb 504287
 Hand Tight- Hose Barb Swivel with 1 1/2" Stem	1/4" Hose Barb 0118S	1/4" Hose Barb 0218S	1/4" Hose Barb 0318S 5/16" Hose Barb 0319S	1/4" Hose Barb 0418S	1/4" Hose Barb 0518S	1/4" Hose Barb 0618S 5/16" Hose Barb 0619S	1/4" Hose Barb 504287S

DISS Male

MEDICAL FITTINGS



Oxygen



Air



Vacuum



Nitrous
Oxide



Nitrogen



WAGD



Carbon Oxide

	Oxygen	Air	Vacuum	Nitrous Oxide	Nitrogen	WAGD	Carbon Oxide
DISS Male to 1/8" NPT Male	0701	0731	0751	0771	0791		
DISS Male to 1/4" NPT Male	0703	0733	0753	0773	0792	0794	8946
 1/4" Hose Barb	0709					0601 Demand Check 0602	
 1/4" Hose Barb with Knob	0710	0740	0760	0780	0795		
 1/4" Hose Barb with Knob & Check	0711	0741	0761	0781	0796		
 1/4" Hose Barb with Knob Swivel	0710S	0740S	0760S	0780S	0795S		
 1/4" Hose Barb with Knob & Check Swivel	0711S	0741S	0761S	0781S	0796S		

DISS Male

Demand Check Valve



Oxygen



Air



Vacuum



Nitrous
Oxide



Nitrogen



WAGD



Carbon Oxide

DISS Male with Demand Check Valve to 1/8" NPT Male	0712	0742		0782	0797		
DISS Male with Demand Check Valve to 1/4" NPT Male	0713	0743	0763	0783	0798	0799	8947

One Way Check Valve



O2, DISS Male x 1/4" NPT Male

508373



Air, DISS Male x 1/4" NPT Male

508374



1/4" NPT Female x 1/4" NPT Male

508375



Oxygen DISS 90° Elbow
Male x Hex Nut

1" Stem
0705
1 1/2" Stem
0706



Oxygen DISS Male x
DISS Male

0704



Oxygen DISS Male x
1/8" NPT Female

0702



Oxygen Power Take-Off
(1/8" NPT)

0150



Air Power Take-Off
(1/8" NPT)

0151

Puritan-Bennett Style Quick Connect



Oxygen



Air



Vacuum



Nitrous Oxide



WAGD

	4101	4201	4301	4401	4601
1/8" NPT Male					
			4301X		
1/8" NPT Male x 1.5" Long					
	4102	4202	4302	4402	4602
1/8" NPT Female					
	4103	4203	4303	4403	4603
1/4" NPT Male					
	4105				
DISS Male with Knob (O2)					
		4205	4305	4405	4605 (no knob)
DISS Male with Knob					
	4106	4206	4306	4406	4606
DISS Hex Nut					
	4108	4208	4308	4408	4608
DISS Hand Tight					
	1/4" Hose Barb 4111	1/4" Hose Barb 4211	1/4" Hose Barb 4311 5/16" Hose Barb 4312 (no knob)	1/4" Hose Barb 4411	1/4" Hose Barb 4611 (no knob) 5/16" Hose Barb 4615 (no knob)
Hose Barb with Knob					
	1/4" Hose Barb 4111S	1/4" Hose Barb 4211S	1/4" Hose Barb 4311S 5/16" Hose Barb 4312S	1/4" Hose Barb 4411S	1/4" Hose Barb 4611S 5/16" Hose Barb 4615S
Hose Barb Swivel					
	4114				
DISS Male with Knob and Check Valve (O2)					
		4214	4314	4414	4614 (no knob)
DISS Male with Knob and Check Valve					

MEDICAL FITTINGS

Puritan-Bennett Style Check Unit (Female Coupler)




Oxygen



Air



Vacuum

	Oxygen	Air	Vacuum
 1/8" NPT Female	4150	4250	4350
 1/8" NPT Male	4151	4251	4351
 1/4" NPT Male	4153	4253	4353
 Hose Barb	1/4" Hose Barb 4154	1/4" Hose Barb 4254	1/4" Hose Barb 4354 5/16" Hose Barb 4357
 DISS Male (O2)	4155		
 DISS Male		4255	4355
 DISS Hex Nut	4156	4256	4356
 DISS Hand Tight	4158	4258	4358

Schrader Style Quick Connect

MEDICAL FITTINGS



Oxygen



Air



Vacuum



Nitrous Oxide

	6101	6201	6301	6401
1/8" NPT Male				
			6301X	
1/8" NPT Male x 1.5" Long				
	6102	6202	6302	6402
1/8" NPT Female				
	6103	6203	6303	6403
1/4" NPT Male				
	6105			
DISS Male (O2)				
		6205	6305	6405
DISS Male				
	6106	6206	6306	6406
DISS Hex Nut				
	6108	6208	6308	6408
DISS Hand Tight				
	1/4" Hose Barb 6111	1/4" Hose Barb 6211	1/4" Hose Barb 6311	1/4" Hose Barb 6411
Hose Barb				
	6114			
DISS Male and Check Valve (O2)				
		6214	6314	6414
DISS Male and Check Valve				

Schrader Piston Identification



Hex Piston
(standard)



Swivel Piston
(add "S" after part number)

Schrader Style Check Unit (Female Coupler)



Oxygen



Air






Vacuum



Nitrous Oxide



Air/Nitrogen

	Oxygen	Air	Vacuum	Nitrous Oxide	Air/Nitrogen
 1/8" NPT Female	6150	6250	6350	6450	
 1/8" NPT Male	6151	6251	6351	6451	
 1/4" NPT Male	6153	6253	6353	6453	
 Hose Barb	1/4" Hose Barb 6154	1/4" Hose Barb 6254	1/4" Hose Barb 6354	1/4" Hose Barb 6454	
 DISS Male (O2)	6155				
 DISS Male		6255	6355	6455	
 DISS Hex Nut	6156	6256	6356	6456	Schrader Air Coupler x Nitrogen DISS 6556
 DISS Hand Tight	6158	6258	6358	6458	Schrader Air Coupler x Nitrogen DISS 6558

Medstar Style Quick Connect



Oxygen



Air



Vacuum

MEDICAL FITTINGS

	8101	8201	8301
1/8" NPT Male			
			8301X
1/8" NPT Male x 1.5" Long			
	8102	8202	8302
1/8" NPT Female			
	8103	8203	8303
1/4" NPT Male			
	8105		
DISS Male (O2)			
		8205	8305
DISS Male			
	8106	8206	8306
DISS Hex Nut			
	8108	8208	8308
DISS Hand Tight			
	1/4" Hose Barb 8111	1/4" Hose Barb 8211	1/4" Hose Barb 8311
Hose Barb			
	8111S	8211S	8311S
Hose Barb Swivel			
	8114MS		
DISS Male with Check Valve (O2)			
		8214	8314
DISS Male with Check Valve			

Oxequip Style Quick Connect

MEDICAL FITTINGS



Oxygen



Air



Vacuum

	Oxygen	Air	Vacuum
	5101	5201	5301
			5301X
	5102	5202	5302
	5103	5203	5303
	5105		
		5205	5305
	5106	5206	5306
	5108	5208	5308
	1/4" Hose Barb 5111	1/4" Hose Barb 5211	1/4" Hose Barb 5311
	5114		
		5214	5314

Y Block Assemblies and DISS Y Assemblies

MEDICAL FITTINGS

	OXYGEN	AIR	VACUUM	NITROUS OXIDE	NITROGEN
 Ohmeda Style	2152Y	2252Y	2352Y	2452Y	
with Winged Quick Connect	2152WY	2252WY	2352WY	2452WY	
 Compact Ohmeda Style	2162Y	2262Y	2362Y	2462Y	
with Winged Quick Connect	2162WY	2262WY	2362WY	2462WY	
 Chemetron Style	3152Y	3252Y	3352Y	3452Y	
 Puritan-Bennett Style	4152Y	4252Y	4352Y		
 Schrader Style	6152Y	6252Y	6352Y	6452Y	
 DISS Y Assemblies	7211	7212	7213	7214	7240

T Branch Assemblies and DISS T Assemblies

OXYGEN

AIR

VACUUM

NITROUS
OXIDE

 Ohmeda Style	2152T	2252T	2352T	2452T
with Winged Quick Connect	2152WT	2252WT	2352WT	2452WT
 Ohmeda Style Offset Right	2152T-R	2252T-R	2352T-R	2452T-R
with Winged Quick Connect	2152WT-R	2252WT-R	2352WT-R	2452WT-R
 Ohmeda Style Offset Left	2152T-L	2252T-L	2352T-L	2452T-L
with Winged Quick Connect	2152WT-L	2252WT-L	2352WT-L	2452WT-L
 Compact Ohmeda Style	2162T	2262T	2362T	2462T
with Winged Quick Connect	2162WT	2262WT	2362WT	2462WT
 Compact Ohmeda Style Offset Right	2162T-R	2262T-R	2362T-R	2462T-R
with Winged Quick Connect	2162WT-R	2262WT-R	2362WT-R	2462WT-R
 Compact Ohmeda Style Offset Left	2162T-L	2262T-L	2362T-L	2462T-L
with Winged Quick Connect	2162WT-L	2262WT-L	2362WT-L	2462WT-L

T Branch Assemblies and DISS T Assemblies

MEDICAL FITTINGS

OXYGEN

AIR

VACUUM

NITROUS
OXIDE



Chemetron Style

3152T

3252T

3352T

3452T



Chemetron Style
Offset Right

3152T-R

3252T-R

3352T-R

3452T-R



Chemetron Style
Offset Left

3152T-L

3252T-L

3352T-L

3452T-L



Puritan-Bennett Style

4152T

4252T

4352T



Puritan-Bennett Style
Offset Right

4152T-R

4252T-R

4352T-R



Puritan-Bennett Style
Offset Left

4152T-L

4252T-L

4352T-L

T Branch Assemblies and DISS T Assemblies

	OXYGEN	AIR	VACUUM	NITROUS OXIDE
 Schrader Style	6152T	6252T	6352T	6452T
 Schrader Style Offset Right	6152T-R	6252T-R	6352T-R	6452T-R
 Schrader Style Offset Left	6152T-L	6252T-L	6352T-L	6452T-L
 DISS T Assemblies	7111T	7112T		7114T
 DISS T Assemblies Offset Right	7111T-R	7112T-R		7114T-R
 DISS T Assemblies Offset Left	7111T-L	7112T-L		7114T-L

T Branch, Y Block, Hex/Quad Block and Rail Blocks

MEDICAL FITTINGS



T Branch



Y Block



Hex Block (3 outlet)




Quad Block (4 outlet)

1/8" NPT Standard
7100
1/8" NPT Offset
7100L/R

1/8" NPT Female
7201
1/4" NPT Female
7200

1/8" NPT Female
9054
1/4" NPT Female
9052

1/4" NPT Female
9051



Hex Block
(1/8" NPT Female Inlet)

7217

7218

7219



Quad Block Assembly
(1/8" NPT Female Inlet)

7220

7221


7222

Rail Block Assemblies
All Rail Block Assemblies include Wall Bracket and Plate

7223

7224

7225



Rail Block Assembly
(DISS Male Inlet)

Replacement Parts
7226 — Mounting Bracket
9056 — Mounting Plate
9054 — Hex Block, 3 outlets
9058 — Screw (2 required)

Extension Manifolds, Pole Clamp and Ferrules

OXYGEN

AIR

VACUUM

Extension Manifold

Features:

Economical expansion of single wall outlet.

Permanently mounts to wall.

2 foot conductive hose with 1/8" NPT Female inlet is standard.

Outlets are DISS Male with Check Valves.

	7230	7233	
Extension Manifold-2 Outlet			
			7236
Extension Manifold-2 Outlet with Swivel			
	7231	7234	
Extension Manifold-3 Outlet			
	7232	7235	
Extension Manifold-4 Outlet			
	505568		
Oxygen Manifold			

Pole Clamp

Features:

1/8" NPT Female holes (top and front) - Comes with 1 plug on bottom

	7227A shown	1" Pole 7216
	7227B	1-1/2" Pole 7227
Pole Clamp		
	1/4" Hose	5/16" Hose

Ferrule

	3/4" Long	1" Long	3/4" Long	1" Long
	9030	507712	503780	507501
Ferrule				
	Add -25 or -50 at end of Part Number for those quantities			

Tubing Nipples

Method of Cleaning

Part Number	Autoclave	Ethylene Oxide	Chlorine Dioxide	Hydrogen Peroxide	Gamma Ray	Electron Beam	Vaporized Hydrogen Peroxide (VHP)	Quaternary Cleaners	Sodium Hypochlorite "Bleach Solutions"
0154/0154C		X	X	X	X	X	X	X	X
0158/0158C		X	X	X	X	X	X	X	X
0159/0159C	X	X	X	X	X	X	X	X	X
0159K/0159KC	X	X	X	X	X	X	X	X	X
0115	X	X	X	X	X	X	X	X	X

* Please Remove O-Ring for Autoclaving. Clean all Tubing Nipples after each patient use. *

 0154 0154C Hand Tight Tubing Nipple	Oxygen	0154
	Oxygen with Chain	0154C
 0158 0158C Wing Nut Tubing Nipple	Oxygen	0158
	Oxygen with Chain	0158C
 0159 0159C Hex Nut Tubing Nipple	Oxygen	0159
	Oxygen with chain	0159C
 0159K 0159KC Knurled Nut Tubing Nipple	Oxygen	0159K
	Oxygen with chain	0159KC
 Disposable Tubing Nipple	Clear	0119
	Yellow	0120
	Black	0121
	White	0122
	Green	0124
*Note: sold in packs of 25 and 50. Add -25 or -50 at end of part number for those quantities		
 9073 9090 Tubing Nipple (for Vacuum Regulator)	Straight	9073
	90°	9090
 Chain Ring Assembly	6" Chain with Closed Ring	8004
	6" Chain with Open Ring	8005

MEDICAL FITTINGS

Medical Grade Hose

Features:

Kink resistant and light weight.

FDA medical grade conductive hose.

Reinforced with tough nylon braiding.

Conveniently packaged in 250 foot rolls or as noted.

Color coded to conform to CGA-C9 standards.

Maximum working pressure to 200 psi @ 70°F.

Meets NFPA 99 requirements for conductivity (conductive hose only).

Keep out of sunlight and in a climate controlled atmosphere

Non-Conductive Hose

Oxygen	Green		(1/4" ID-250' Roll)	9020-250
Air	Yellow		(1/4" ID-250' Roll)	9021-250
Vacuum	White		(1/4" ID-250' Roll)	9022-250
Nitrous Oxide	Blue		(1/4" ID-250' Roll)	9023-250
Nitrogen	Black		(1/4" ID-250' Roll)	9028-250
Carbon Dioxide	Gray		(1/4" ID-250' Roll)	9033-250
WAGD	Purple		(5/16" ID-250' Roll)	9038-250
	Clear		(1/4" ID-100' Roll)	504202-100

Conductive Hose

Oxygen	Green		(1/4" ID-250' Roll)	9120-250
Air	Yellow		(1/4" ID-250' Roll)	9121-250
Vacuum	White		(1/4" ID-250' Roll) (5/16" ID-250' Roll)	9122-250 9135-250
Nitrous Oxide	Blue		(1/4" ID-250' Roll)	9123-250
Nitrogen	Black		(1/4" ID-250' Roll)	9128-250
Carbon Dioxide	Gray		(1/4" ID-250' Roll)	9125-250
WAGD	Purple		(1/4" ID-250' Roll)	9124-250
ISO Air	Black & White Striped		(1/4" ID-250' Roll)	9126-250

Custom hose assembly configurator is available at: <https://precisionmedical.com/hose-assemblies>



Crimping Tool
9130



Hose Retractor

1/4" Single Clamp	HR-1
1/4" Double Clamp	HR-2
5/16" Single Clamp	HR-3
5/16" Double Clamp	HR-4

Hose Tethers

MEDICAL FITTINGS

OXYGEN

AIR

VACUUM

NITROUS
OXIDE

WAGD



Ohmeda Style Quick Connect

2105HT

2205HT

2305HT

2405HT

2605HT



Chemetron Style

3105HT

3205HT

3305HT

3405HT

3605HT



Puritan-Bennett Style

4105HT

4205HT

4305HT

4405HT

4605HT



Oxequip Style

5105HT

5205HT

5305HT



Schrader Style

6105HT

6205HT

6305HT

6405HT



Medstar Style

8105HT

8205HT

8305HT

NPT Pipe Fittings

	1/8" NPT	1/4" NPT
 Pipe Tee	7030	7031
 Street Tee	7032	7033
 Female Tee	7034	7035
 Plug	7019	7017
 Drop Ear Tee 7081 Elbow 7082		
 Elbow	7092	

	1/8" NPT	1/4" NPT
 Plug	7013	7012
 45° Street Elbow	7056	7057
 90° Elbow	7058	7059
 90° Street Elbow	7063	7065
 Master Flow Valve	7064	7045
 Cross	7036	7037

NPT Pipe Fittings

MEDICAL FITTINGS

 Hex Nipple	1/8" NPT Male x 1/8" NPT Male	7020
	1/4" NPT Male x 1/4" NPT Male	7021
	1/8" NPT Male x 1/4" NPT Male	7022
 Connector	1/8" NPT Female x 1/8" NPT Male	1017
	1/8" NPT Female x 1/4" NPT Male	7024
	1/4" NPT Female x 1/8" NPT Male	7025
	1/4" NPT Female x 1/4" NPT Male	7026
	1/8" NPT Female x 1/8" NPT Male x 2.5" Long	8955
 Coupler	1/8" NPT Female x 1/8" NPT Female	7027
	1/4" NPT Female x 1/4" NPT Female	7028
	1/4" NPT Female x 1/8" NPT Female	7029
 Hose Barb Connector	1/8" NPT Male x 1/8" Diameter	7038
	1/8" NPT Male x 3/16" Diameter	7039
	1/8" NPT Male x 1/4" Diameter	7040
	1/4" NPT Male x 1/4" Diameter	7043
	1/8" NPT Male x 5/16" Diameter	8952
	1/8" NPT Female x 1/4" Diameter	7046
	1/4" NPT Female x 1/4" Diameter	7049
	1/4" NPT Female x 5/16" Diameter	8951
 Swiveling Hose Barbs	1/8" NPT Male x 5/16" Diameter	507765
	1/8" NPT Male x 1/4" Diameter	507764
 Pipe Nipple	1/8" NPT Male x 0.75" Long	7044
	1/8" NPT Male x 1.5" Long	7050
	1/8" NPT Male x 2" Long	7051
	1/8" NPT Male x 2.5" Long	7052
	1/4" NPT Male x 1.5" Long	7053
	1/4" NPT Male x 2" Long	7054
	1/4" NPT Male x 2.5" Long	7055
	1/8" NPT Male x 3" Long	7062
 Cap	1/8" NPT Female	7018
 Hose Splicer	1/4" Inside Diameter	7041

Neb Drives and Rescue Flows

MEDICAL FITTINGS

Neb Drives

		OXYGEN AIR	
 Ohmeda	6 to 7 l/min	8126	8226
	7 to 8 l/min	8127	8227
	8 to 9 l/min	8128	8228
	10 to 11 l/min	8129	8229
 Chemetron	6 to 7 l/min	8136	8236
	7 to 8 l/min	8137	8237
	8 to 9 l/min	8138	8238
	10 to 11 l/min	8139	8239
 DISS Hand Tight	6 to 7 l/min	8116	8216
	7 to 8 l/min	8117	8217
	8 to 9 l/min	8118	8218
	10 to 11 l/min	8119	8219
 Puritan-Bennett with Knob	6 to 7 l/min	8146	8246
	7 to 8 l/min	8147	8247
	8 to 9 l/min	8148	8248
	10 to 11 l/min	8149	8249
 Schrader	6 to 7 l/min	8166	8266
	7 to 8 l/min	8167	8267
	8 to 9 l/min	8168	8268
	10 to 11 l/min	8169	8269
 Oxequip with Knob	6 to 7 l/min	8156	8256
	7 to 8 l/min	8157	8257
	8 to 9 l/min	8158	8258
	10 to 11 l/min	8159	8259
 Medstar	6 to 7 l/min	8176	8276
	7 to 8 l/min	8177	8277
	8 to 9 l/min	8178	8278
	10 to 11 l/min	8179	8279

Rescue Flows

 OXYGEN Oxygen - 14 to 16 l/min	
DISS Hand Tight	8115
Ohmeda	8125
Chemetron	8135
Puritan-Bennett	8145
Oxequip	8155
Schrader	8165
Medstar	8175

 OXYGEN Oxygen - 24 to 26 l/min	
DISS Hand Tight	8114
Ohmeda	8124
Chemetron	8134
Puritan-Bennett	8144
Oxequip	8154
Schrader	8164
Medstar	8174

Condensation Traps

	Inlet DISS Female Oxygen x Outlet DISS Male Oxygen	PM15-4
	Inlet DISS Female Air x Outlet DISS Male Air	PM15-41
	Inlet DISS Female Air x Outlet DISS Male Air 90° Down	PM15-43
	Inlet DISS Female Oxygen 90° Up x Out- let DISS Male Oxygen 90° Down	PM15-44
	Inlet DISS Male Air 90° Down x Outlet DISS Female Air 90° Up	PM15-45

Replacement Parts

	507588
	507587
	1271

Ways to place an order:



Contact your Precision Medical sales representative



Visit precisionmedical.com



Contact Customer Service:
orders@precisionmedical.com
1-800-272-7285

Do you need a special configuration?

Custom configurations, international specifications, special ratings and color options are available.

Contact Customer Service at 1-800-272-7285 for custom requests and more information.